



The *Empress* Safety Huber Needle from Osiris Medical is designed with both the patient and clinician in mind. The Empress maximizes patient comfort and facilitates clinician use and safety. The patented design requires no change in technique and it's a snap to use!

### ***ENHANCING PATIENT COMFORT . . .***

***“Small and light-weight”***, patient forgets it's there!

***“Flexibly designed”*** eliminating patient wear issues and discomfort inherent with ridged type devices.

***“Lowest profile”*** device in market no unsightly bulges under clothing. The barrel lock keeps device in place preventing tenting.

***“Integrated-like comfort pad”*** for comfortable resting of device on patient skin.

***“Siliconized needle”*** facilitating insertion and reducing patient discomfort associated with needle insertion and removal.

Patient trauma reduced through ***“Controlled removal”*** from port.

### ***FACILITATING CLINICIAN USE AND SAFETY . . .***

***“Device is Intuitive”*** requiring no significant training for the experienced clinician.

***“Siliconized needle”*** supporting clinician insertion and removal.

***“Controlled removal”*** eliminates needle “bounce-back” and inadvertent needle sticks to both clinician and patient.

Device allows clinician to reposition, if necessary, without unplanned lock activation.

***“Patient Portholes”*** allows clinician to visually observe injection site for signs of patient redness and irritation.

Visual and audible (snap) confirms safety activation.